


Upholstery Colour Options

We use Studio Encore upholstery made by Austex Commercial Fabrics on all our products, because it offers long lasting quality in a wide range of colours. The vinyl fabric is flame resistant, UV resistant and is treated with a Protech coating, which helps to guard against stains and inhibits the growth of bacteria in healthcare environments. It has also been formulated to withstand the hardening and breakdown of the PVC surface due to perspiration and oils.

We offer the standard colour range with all our products. The premium colour range is available on request for an additional price.

Standard colour range



Premium colour range



*Not all colours shown. More can be seen at www.austexindustrialfabrics.com.au.

The colours shown here may not be a true representation of the actual product and must be viewed as a guide only.

Quality Assurance

Lead and chromium free

Flame Resistant properties: Tested to Australian standards (AS 1530.3)

Spread of flame: 0

I.M.O Standards A. 652 (16)

F.A.S. - 25.853



Australian
Medical Couches

Unit 2, 2 Network Drive, Carrum Downs VIC 3201 AUSTRALIA
P: +61 3 9708 2661 F: +61 3 9708 2617
E: info@austmc.com austmc.com