




Forme **Medical**

evolving with you



PRODUCT CATALOGUE



OUR PRODUCT RANGE

We value being innovative, adaptable and accommodating in the development of our products and our customer service.

Our current product range has subsequently expanded to evolve with you. Following are patient positioning equipment options to service the following areas:

- GENERAL EXAMINATION
- ULTRASOUND
- TREATMENT ROOMS
- GYNAECOLOGY & OBSTETRICS
- VASCULAR
- CARDIOLOGY
- DERMATOLOGY
- ENT
- PODIATRY
- COSMETIC SURGERY

Tel. 03 9708 2661
www.formemedical.com.au

ARGON

AMC 1020



- ECONOMY BASED 2 SECTION BED
- ELECTRIC HEIGHT & BACKREST ADJUSTMENT
- FOUR LOCKING BRAKE CASTORS
- HAND CONTROL



ONYX

AMC 2510



- 2 SECTIONS
- ELECTRIC HEIGHT & BACKREST ADJUSTMENT
- FOUR LOCKING BRAKE CASTORS
- HAND CONTROL



TANZANITE

AMC 2520



60 DEGREE TILT TABLE

ELECTRIC TILT, HEIGHT & BACKREST ADJUSTMENT

CENTRAL LOCKING

BATTERY BACK UP



MALACHITE

AMC 2540



ECHO CARDIO DROP SECTION

ELECTRIC HEIGHT & BACKREST ADJUSTMENT

FOUR LOCKING BRAKE CASTORS

FOOT CONTROL



QUARTZ

AMC 2530



CUT OUT SECTION FOR IMPROVED PATIENT ACCESS

ELECTRIC HEIGHT & BACKREST ADJUSTMENT

FOUR LOCKING BRAKE CASTORS

FOOT CONTROL



OPAL

AMC 2550



3 SECTIONS

ELECTRIC HEIGHT, LEGREST & BACKREST ADJUSTMENT

FOUR LOCKING BRAKE CASTORS

HAND CONTROL



CITRINE

AMC 2570



2 SECTIONS WITH 12 DEGREE TRENDLENBURG TILT
ELECTRIC TILT, HEIGHT & BACKREST ADJUSTMENT
FOUR LOCKING BRAKE CASTORS
HAND CONTROL



TOPAZ

AMC 2120



4 SECTION CHAIR WITH INDEPENDENT FOOT SECTIONS
ELECTRIC HEIGHT & BACKREST ADJUSTMENT
FOUR LOCKING BRAKE CASTORS
FOOT CONTROL



SPHENE

AMC 2110



3 SECTION SPECIALISTS CHAIR WITH SHAPED BACK
ELECTRIC HEIGHT, LEGREST & BACKREST ADJUSTMENT
FOUR LOCKING BRAKE CASTORS
HAND CONTROL



SAPPHIRE

AMC 2130



FULLY CUSTOMISABLE PROCEDURE CHAIR / BED
ELECTRIC HEIGHT & BACKREST ADJUSTMENT
FOUR LOCKING BRAKE CASTORS
FOOT CONTROL



AMETHYST

AMC 2140



3 SECTION CHAIR / BED

ELECTRIC HEIGHT, FOOTREST & BACKREST ADJUSTMENT

FOUR LOCKING BRAKE CASTORS

HAND CONTROL



OPERATOR STOOL

AMC OP



HEIGHT ADJUSTABLE

FIVE CASTORS

RANGE OF COLOURS



ERGO CHAIR

AMC ER

HEIGHT ADJUSTABLE

HEIGHT & TILT ADJUSTABLE BACKREST

RANGE OF COLOURS



WHY FORME MEDICAL?

If you are in the market for a new treatment chair, table or examination couch you want to be confident that the product you buy has the functionality and long lasting quality that you need. Being a capital equipment item that you or your staff members will be using all day every day, you need to be sure that the equipment improves the quality of life for both the staff and the patient.

Tel. 03 9708 2661

www.formemedical.com.au

PURA TREATMENT CHAIR

DIA AX



- MULTIFUNCTIONAL REMOVABLE ARMRESTS
- RECLINES TO TRENDELENBURG & SUPINE POSITION
- INDIVIDUAL LOCK BRAKING
- PATIENT ADJUSTING FOOT SUPPORT



- NON-HEIGHT ADJUSTABLE
- 3 MOTORS CONTROLLING BACK, SEAT & FOOT SECTIONS
- BATTERY BACK-UP
- EMERGENCY BUTTON POSITIONING



PURA TREATMENT CHAIR

DIA BX



- MULTIFUNCTIONAL REMOVABLE ARMRESTS
- RECLINES TO TRENDELENBURG & SUPINE POSITION
- INDIVIDUAL LOCK BRAKING
- ADVANCED PATIENT ADJUSTING FOOT SUPPORT



- HEIGHT ADJUSTABLE
- 4 MOTORS CONTROLLING LIFT, BACK, SEAT & FOOT SECTIONS
- BATTERY BACK-UP
- EMERGENCY BUTTON POSITIONING



PURA TREATMENT CHAIR

DIA CX



MULTIFUNCTIONAL REMOVABLE ARMRESTS
RECLINES TO TRENDELENBURG & SUPINE POSITION
CENTRAL LOCK BRAKING
ADVANCED PATIENT ADJUSTING FOOT SUPPORT

HEIGHT ADJUSTABLE
4 MOTORS CONTROLLING LIFT, BACK, SEAT & FOOT SECTIONS
BATTERY BACK-UP
EMERGENCY BUTTON POSITIONING



ONE DAY SURGERY CHAIR

DIA FX










MULTIFUNCTIONAL ARMRESTS THAT ACTS AS A
SIDE RAIL & DRINK HOLDER
EMERGENCY FOOT RELEASE FOR HANDSFREE OPERATION
CENTRAL LOCK BRAKING, WITH DIRECTIONAL LOCK
ADVANCED PATIENT ADJUSTING FOOT SUPPORT

HEIGHT ADJUSTABLE
4 MOTORS CONTROLLING LIFT, BACK, SEAT & FOOT SECTIONS
BATTERY BACK-UP
EMERGENCY BUTTON POSITIONING
- TRENDELENBURG & SUPINE POSITION



TREATMENT CHAIRS

We have sourced the globe to bring the very best in medical treatment chair technology to our local market. Introducing the Pura range.

| | |
|---|-----------------------------|
|  | RENAL DIALYSIS |
|  | CHEMOTHERAPY / DAY ONCOLOGY |
|  | BLOOD COLLECTION |
|  | DAY RECOVERY |
|  | ICU |
|  | DAY SURGERY |
|  | SPECIALIST APPLICATIONS |

SPECIFICATIONS

| | WIDTH mm | LENGTH mm | BACK REST TILT Degrees | HEIGHT ADJUSTABLE mm | LIFTING CAPACITY kgs | AUSTRALIAN MADE | WARRANTY | |
|-----------|-------------|--------------|------------------------------|----------------------------|----------------------------|--------------------|----------|------------|
| | | | | | | | Frame | Electrical |
| ARGON | 600 | 1950 | 75 | 550-1060 | 200 | ✗ | 5 years | 2 years |
| ONYX | 610 | 1970 | 75 | 550-1000 | 200 | ✓ | Lifetime | 5 years |
| TANZANITE | 610 | 2000 | 70 | 550-1000 | 200 | ✓ | Lifetime | 5 years |
| QUARTZ | 700 | 1970 | 75 | 550-1000 | 200 | ✓ | Lifetime | 5 years |
| MALACHITE | 670 | 1970 | 75 | 550-1000 | 200 | ✓ | Lifetime | 5 years |
| OPAL | 610 | 1970 | 75 | 550-1000 | 200 | ✓ | Lifetime | 5 years |
| CITRINE | 670 | 1970 | 75 | 550-1000 | 200 | ✓ | Lifetime | 5 years |
| SPHENE | 610 | 1850 | 75 | 450-1000 | 200 | ✓ | Lifetime | 5 years |
| TOPAZ | 610 | 1870 | 75 | 550-1000 | 200 | ✓ | Lifetime | 5 years |
| SAPPHIRE | 610 | 1350 | 75 | 450-1000 | 200 | ✓ | Lifetime | 5 years |
| AMETHYST | 610 | 1850 | 75 | 450-1000 | 200 | ✓ | Lifetime | 5 years |

| | DIAMETER mm | HEIGHT ADJ. mm | LIFTING CAPACITY kgs | AUSTRALIAN MADE | WARRANTY |
|----------------|----------------|----------------------|----------------------------|--------------------|----------|
| OPERATOR STOOL | 350 | 490-690 | 120 | ✓ | 1 year |
| ERGO STOOL | 350 | 490-690 | 120 | ✓ | 1 year |

| | WIDTH OF REST AREA mm | LENGTH OF REST AREA mm | TTL WIDTH W/OUT ARMRESTS mm | TTL WIDTH WITH ARMRESTS mm | HEIGHT ADJUSTMENT mm | LIFTING CAPACITY kgs |
|-------------------|-----------------------------|------------------------------|-----------------------------------|----------------------------------|----------------------------|----------------------------|
| PURA DIA AX** | 590 | 2100 | 680 | 910 | 530 | 180 |
| PURA DIA BX** | 590 | 2100 | 680 | 910 | 515-765 | 180 |
| PURA DIA CX** | 590 | 2100 | 730 | 910 | 515-765 | 200 |
| ONE DAY SURGERY** | 590 | 2100 | 680 | 730 | 550-800 | 200 |

**Made in Czech Republic * 2 Year Warranty



11 Frankston Gardens Drive
Carrum Downs 3201, VIC AUS
Tel. (03) 9708 2661 Fax. (03) 9708 2617
Email. info@formemedical.com.au www.formemedical.com.au